


## Unambiguous Molecular Identification and Quantification

**Molecular Rotational Resonance (MRR)** simultaneously provides unambiguous identification and quantification of small molecules by measuring their unique spectral fingerprint.

MRR's chemical specificity easily discriminates between all types of isomers including regioisomers (shown in figure below), diastereomers, enantiomers (providing EE and AC), and isotopic species with site-specific information.

## Structural Validation

MRR identifies molecules through the characterization of their rotational spectra, where subtle differences in structure produce significant changes in their spectra. Quantum chemistry calculations allow for the unambiguous confirmation of the molecular identity, without a pure reference standard.



## Routine Analysis

High sample throughput is achieved by the speed of MRR analysis with minimal sample preparation.

Additionally the high resolution, sensitivity, and specificity of MRR eliminates the need for front-end separation allowing for the direct analysis of complex mixtures, enabling online reaction monitoring.

Quantitation is achieved without reference standards as the signal is directly proportional to the sample concentration.

### Several recent publications provide deeper insight into the utility of MRR:

**Chiral Analysis of Linalool, an Important Natural Fragrance and Flavor Compound, by Molecular Rotational Resonance Spectroscopy**

[BrightSpec \(Sonstrom, Cannon, Neil\): \*Symmetry\*, 2022, 14, 917](#) [OPEN ACCESS]

**Chiral Analysis of Pantolactone with Molecular Rotational Resonance Spectroscopy**

[BASF, Pate and co-workers: \*Chirality\*, 2022, 34, 114-135](#) [OPEN ACCESS]

**Direct regioisomer analysis of crude reaction mixtures via molecular rotational resonance (MRR) spectroscopy**

[Merck, Pate and co-workers: \*Chem. Sci.\*, 2020, 11, 6332-6338](#)

**Rapid quantification of isomeric and dehalogenated impurities in pharmaceutical raw materials using MRR spectroscopy**

[GSK, Pate and co-workers: \*Journal of Pharmaceutical and Biomedical Analysis\*, 2020, 10, 113474](#)

## BrightSpec Product Portfolio

BrightSpec offers a complete line of MRR spectrometers. The Broadband MRR delivers comprehensive sample profiling and characterization of molecular structures, including new compounds, without a reference standard.

The isoMRR provides routine analysis, screening and quantitation, and direct chiral or achiral process monitoring.

Additionally, BrightSpec offers a line of MRR spectrometers designed for Educational and Fundamental research. Please reach out for more information on what BrightSpec has to offer.



### Broadband MRR Structural Elucidation and Research

<b>Typical Applications</b>	Structural Elucidation Enantiomeric Excess Quantitative Mixture Analysis
<b>Frequency Range</b>	2–8 GHz, 6–18 GHz
<b>Analyte Range</b>	100–400 amu
<b>Analysis Time</b>	60+ minutes
<b>Sample Type</b>	Solids, Liquids and Gases
<b>Sample Volume</b>	5–50 mg (typical)

\*Specifications are sample dependent

### isoMRR Routine and Targeted Analysis

<b>Typical Applications</b>	Isomeric Impurity Quantification Enantiomeric Excess Continuous Process Monitoring
<b>Frequency Range</b>	6–18 GHz
<b>Analyte Range</b>	100–400 amu
<b>Analysis Time</b>	5–20 minutes
<b>Sample Type</b>	Solids, Liquids and Gases
<b>Sample Volume</b>	0.5–5 mg (typical)

\*Specifications are sample dependent

[sales@brightspec.com](mailto:sales@brightspec.com)

[www.brightspec.com](http://www.brightspec.com)

